

TENS 3000™

The TENS 3000™ analog unit uses a microprocessor, which gives exact dial performance. Benefits includes maximum full strength power of 80 mA. The TENS 3000™ is as accurate as a digital unit, yet easy to operate. The safety amplitude cap protects the amplitude knobs from accidental bumps or knob movements.



Product Features:

- System includes device, lead wires, four self-adhesive reusable electrodes, 9V battery, hard plastic carrying case and instruction manual
- Power Source: 9-volt battery
- Modes: Standard, Modulate and Burst
- Channels: Dual
- Timer: 30 and 60 minutes or continuous
- Pulse Width: 50-250 Microseconds, Adjustable
- Pulse Rate: 2-150 Hz
- Pulse Amplitude: 0-80 mA each channel
- Waveform: Asymmetrical biphasic square pulse
- One year warranty

Item #	Description	Qty
DT3002	TENS 3000™ Analog Unit, 3 Modes	each