

TENS AA

The TENS AA is a unique digital TENS device that uses AA batteries, has five treatment modes and a large LCD screen. This unit has a recall program and is easy for anyone to use.



Product Features:

- Power Source: two AA batteries
- Modes: Constant, Burst I, Burst II, Modulated Width, Modulated Rate and Width
- Channels: Dual
- Timer: 15, 30, 45, 60 and 90 minutes or continuous
- Pulse Frequency: 20-150 Hz
- Waveform: Symmetrical biphasic or mono phasic
- Two year warranty

Item #	Description	Qty
DA1812	TENS AA, Digital Unit, 5 Modes	each